

Recombinant Human  $\beta$ -Defensin 4,DEFB4

Catalog No: #AP74862



Package Size: #AP74862-1 10ug #AP74862-2 50ug #AP74862-3 500ug #AP74862-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human $\beta$ -Defensin 4,DEFB4
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8WTQ1
Uniprot	Q8WTQ1
GeneID	140596;503618;
Calculated MW	5.99kDa
Target Sequence	EFELDRICGYGTARCRKKCRSQEYRIGRCPNTYACCLRKWDESLNRTKP
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 130mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only