

## Recombinant Human N-Glycosylase,OGG1 (N-GST)

Catalog No: #AP74864



Package Size: #AP74864-1 10ug #AP74864-2 50ug #AP74864-3 500ug #AP74864-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human N-Glycosylase,OGG1 (N-GST)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15527
Uniprot	O15527
GeneID	4968;
Calculated MW	65.07kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSR IAYSKDFETLKVDFLSKLP EMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFK KRIEAIPIQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSDL VPRGSM PARALLPRRMGHRTLASTPALWASIPCRSELRLDLVLP SGQSFRWREQSPAHW SGLADQVWTL TQTEEQ LHCTVYRGDKSQASRPTPDELEAVRKYFQLDVTLAQLYHHWGSVDSHFQEVAQKFQGVRLLRQDPI ECLFSFICSSNNNIARITGMVERLCQAFGPRLIQLDDVYHGFPSLQALAGPEVEAHLRKLGLGYRARYVSASA RAILEEQGGLAWLQQLRESSYEEAHKALCILPGVGTKVADCICLMALDKPQAVPVDVHMWHIAQRDYSWHPT TSQAKGSPQTNKELGNFFRSLWGPYAGWAQAVLFSADLRQSRHAQEPAPAKRRRKGSKGPEG
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 1mM EDTA, 50% Glycerol, pH 7.8.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only