

Recombinant Human Deoxyhypusine Synthase,DHS (C-6His)

Catalog No: #AP74865

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74865-1 10ug #AP74865-2 50ug #AP74865-3 500ug #AP74865-4 1mg

Description

Product Name	Recombinant Human Deoxyhypusine Synthase,DHS (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P49366
Uniprot	P49366
GeneID	1725;
Calculated MW	42.04kDa
Target Sequence	MEGSLEREAPAGALAAVLKHSSTLPPESTQVRGYDFNRGVNYRALLEAFGTTGFQATNFGRAVQQVNAMIEK KLEPLSQDEDEDQHADLTQSRRLTSTIFLGYTSNLISSGIRETIRYLVQHNMVDVLVTTAGGVEEDLIKCLAPTY LGEFSLRGKELRENGINRIGNLLVPNENYCKFEDWLMPILDQMVMEQNTEGVKWTPSKMIARLGKEINNPESV YYWAQKNHIPVFSPALTDGSLGDMIFFHSYKNPGLVLDIVEDLRLINTQAIFAKCTGMILGGGVVKKHHIANANL MRNGADYAVYINTAQEFDGSDSGARPDEAVSWGKIRVDAQPVKVYADASLVFPLLVAETFAQKMDAFMHEK NEDLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 20% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only