

Recombinant E. coli Ecotin (C-6His)

Catalog No: #AP74866



Package Size: #AP74866-1 10ug #AP74866-2 50ug #AP74866-3 500ug #AP74866-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E. coli Ecotin (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	E.coli
Accession No.	P23827
Uniprot	P23827
GeneID	946700;
Calculated MW	19.2kDa
Target Sequence	MKTILPAVLFAAFATTSAAWAAESVQPLEKIAPYPQAEKGMKRQVIQLTPQEDESTLKVLLIGQTLEVDCNLHRL GGKLENKTLEGWGYDYYVFDKVSSPVSTMMACPDGKKEKKFVTAYLGDAGMLRYNSKLPVVYTPDNVDVK YRVWKAEEKIDNAVVRLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 300mM NaCl, pH 8.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only