

# Recombinant Human Fatty Acid-Binding Protein 2, FABP2, I-FABP ((N, C-6His)

Catalog No: #AP74870

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74870-1 10ug #AP74870-2 50ug #AP74870-3 500ug #AP74870-4 1mg

## Description

Product Name	Recombinant Human Fatty Acid-Binding Protein 2, FABP2, I-FABP ((N, C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P12104
Uniprot	P12104
GeneID	2169;
Calculated MW	18.44kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAFDSTWKVDRSENYDKFMEKMGVNIVKRKLAHDNLKLTITQEGNKFTVK ESSAFRNIEVVFELGVTFNYNLADGTELRGTWSLEGNKLGKFKRTDNGNELNTVREIIGDELVQTYVYEGVEA KRIFKLDLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only