

Recombinant Human Fatty Acid-Binding Protein 5, FABP5, E-FABP (N-6His)

Catalog No: #AP74873

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74873-1 10ug #AP74873-2 50ug #AP74873-3 500ug #AP74873-4 1mg

Description

Product Name	Recombinant Human Fatty Acid-Binding Protein 5, FABP5, E-FABP (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q01469
Uniprot	Q01469
GeneID	2171;
Calculated MW	17.33kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMATVQQLEGRWRLVDSKGFDEYMKELGVGIALRKMAMAKPDCIITCDGKN LTIKTESTLKTTQFSC TLGKFEETTADGRKTQTVCNFTD GALVQHQEWDGKESTITRKLKDGKLVVECVMMNN VTCTRIYEKVE
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only