

Recombinant Human Fructose-1,6-Bisphosphatase 1,FBPase 1 (E. coli, C-6His)



Catalog No: #AP74875

Orders: order@signalwayantibody.com

Package Size: #AP74875-1 10ug #AP74875-2 50ug #AP74875-3 500ug #AP74875-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fructose-1,6-Bisphosphatase 1,FBPase 1 (E. coli, C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09467
Uniprot	P09467
GeneID	2203;
Calculated MW	37.89kDa
Target Sequence	ADQAPFDTDVNTLTRFVMEEGRKARGTGELTQLLNSLCTAVKAISSAVRKAGIAHLYGIAGSTNVTGDQVKKL DVLNDLVMNMLKSSFATCVLVSEEDKHAIIVEPEKRGKYVVCDFPLDGSSNIDCLVSVGTIFGIYRKKSTDEPS EKDALQPGRNLVAAGYALYGSATMLVLAMDCGVNCFMLDPAIGEFILVDKDKIKKKGIYSLNEGYARDFDPA VTEYIQRKKFPPDNSAPYGARYVGSMAVDVHRTLVIYGGIFLYPANKKSPNGKLRLLYECNPMAYVMEKAGGM ATTGKEAVLDVIPTDIHQRAPVILGSPDDVLEFLKVYKHSQAQVEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, 1mM EDTA, 20% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only