

Recombinant Human Flap Endonuclease 1,FEN-1

Catalog No: #AP74876



Package Size: #AP74876-1 10ug #AP74876-2 50ug #AP74876-3 500ug #AP74876-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Flap Endonuclease 1,FEN-1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P39748
Uniprot	P39748
GeneID	2237;
Calculated MW	42.59kDa
Target Sequence	<p> MGIQGLAKLIADVAPSAIRENDIKSYFGRKVAIDASMSIYQFLIAVRQGGDVLQNEEGETTSHLMGMFYRTIRM MENGIKPVYVFDGKPPQLKSGELAKRSERRAEAEKQLQQAQAAGAEQEVEKFTKRLVKVTKQHNDCKHLLS LMGIPYLDAPSEAEASCAALVKAGKVYAAATEDMDCLTFGSPVLMRHLTASEAKKLPIQEFHLSRILQELGLNQ EQFVDLCILLGSDYCESIRGIGPKRAVDLIQKHKSIEEIVRRDPNKYPVPENWLHKEAHQLFLEPEVLDPESVE LKWSEPNEEELIKFMCGEKQFSEERIRSGVKRLSKSRQGSTQGRLLDDFFKVTGSLSSAKRKEPEPKGSTKKK AKTGAAGKFKRGK </p>
Formulation	Supplied as a 0.2 µm filtered solution of 50 mM Tris, 50mM NaCl, 1mM DDT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only