

## TIP60 (Phospho-Ser86) Polyclonal Antibody

Catalog No: #12237

Package Size: #12237-1 50ul #12237-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TIP60 (Phospho-Ser86) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-TIP60 (S86) Polyclonal Antibody detects endogenous levels of TIP60 protein only when phosphorylated at S86.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human TIP60 around the phosphorylation site of S86.
Target Name	TIP60
Modification	Phospho
Other Names	KAT5; HTATIP; TIP60; Histone acetyltransferase KAT5; 60 kDa Tat-interactive protein; Tip60; Histone acetyltransferase HTATIP; HIV-1 Tat interactive protein; Lysine acetyltransferase 5; cPLA(2)-interacting protein
Accession No.	Swiss-Prot: Q92993NCBI Gene ID: 10524
Uniprot	Q92993
GeneID	10524;
Target Species	human
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

ELISA: 1/40000

Not yet tested in other applications

## Images

Western Blot analysis of Jurkat cells using Phospho-TIP60 (S86) Polyclonal Antibody



Note: This product is for in vitro research use only