

Recombinant Human GABA(A) Receptor-Associated Protein, GABARAP (N-GST)



Catalog No: #AP74878

Orders: order@signalwayantibody.com

Package Size: #AP74878-1 10ug #AP74878-2 50ug #AP74878-3 500ug #AP74878-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GABA(A) Receptor-Associated Protein, GABARAP (N-GST)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q6IAW1
Uniprot	Q6IAW1
Calculated MW	40.21kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSR IAYS KDFETLKVDFLSKLPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIPQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSDL VPRGSMKFVYKEEHPFEKRRSEGEKIRKKYPDRVPVIVEKAPKARIGDLDKKKYLVPSDLTVGQFYFLIRKRIHL RAEDALFFFVNNVIPPTSATMGQLYQEHHEEDFFLYIAYSDES VYGL
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM TrisHCl, 200mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only