

Recombinant Human Grancalcin,GCA (N-GST)

Catalog No: #AP74879



Package Size: #AP74879-1 10ug #AP74879-2 50ug #AP74879-3 500ug #AP74879-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Grancalcin,GCA (N-GST)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P28676
Uniprot	P28676
GeneID	25801;
Calculated MW	50.3kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSRAYSDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIPIQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSDL VPRGSMAYPGYGGGFGNFSIQVPGMQMGQPVPETGPAILLDGYSGPAYSDTYSSAGDSVYTYFSAVAGQDG EVDAAEELQRCLTQSGINGTYSFSLQETCRIMIAMLDRDHTGKMGFNAFKELWAALNAWKENFMTVDQDGSST VEHHELQAIGLMGYRLSPQTLTTIVKRYSKNGRIFDDYVACCVKLRALTDFFRKRDLHQQGSANFIYDDFLQ GTMAI
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 500mM NaCl, pH 8.5 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only