

Recombinant Human Heme Oxygenase 1,HO-1

Catalog No: #AP74881



Package Size: #AP74881-1 10ug #AP74881-2 50ug #AP74881-3 500ug #AP74881-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Heme Oxygenase 1,HO-1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09601
Uniprot	P09601
GeneID	3162;
Calculated MW	29.86kDa
Target Sequence	MERPQPDSMPQDLSEALKEATKEVHTQAENAEFMRNFQKGQVTRDGFKLVMSLYHIYVALEEEIERNKESP VFAPVYFPEELHRKAALEQDLAFWYGPRWQEVIPYTPAMQRYVKRLHEVGRTEPELLVAHAYTRYLGDLSSG QVLKKAQKALDLPSSGEGLAFFTFPNIASATKFKQLYRSRMNSLEMTPAVRQRVIEEAKTAFLLNIQLFEELQE LLTHDTKDQSPSRAPGLRQRASNKVQDSAPVETPRGKPLNT
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only