

Recombinant Human β -Lactamase-Like Protein 2,LACTB2 (N-GST)

Catalog No: #AP74884

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74884-1 10ug #AP74884-2 50ug #AP74884-3 500ug #AP74884-4 1mg

Description

Product Name	Recombinant Human β -Lactamase-Like Protein 2,LACTB2 (N-GST)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q53H82
Uniprot	Q53H82
GeneID	51110;
Calculated MW	59.24kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSRAYSKDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDL VPRGSHMAAVLQRVERLSNRVVRVLGCNPGPMTLQGTNTYLVGTGPRRILIDTGEPAPIEYISCLKQALTEFNT AIQEIVVTHWHRDHSGGIGDICKSINDDTYCIKKLPRNPQREEIIGNGEQYVYLKGDVIKTEGATLRVLYTP GHTDDHMALLEEENAFSGDCILGEGTTFEDLDYDMNSLKELLKIKADIIYPGHGVPVIHNAEAKIQQYISHRNI REQQILTLFRENFEKSFTVMELVKIYKNTPENLHEMAKHNLHLKLEKEGKIFSNTDPDKKWAHL
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only