

Recombinant E. coli Methionine Aminopeptidase, MetAP, MAP (C-6His)

Catalog No: #AP74886

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74886-1 10ug #AP74886-2 50ug #AP74886-3 500ug #AP74886-4 1mg

Description

Product Name	Recombinant E. coli Methionine Aminopeptidase, MetAP, MAP (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	E.coli
Accession No.	P0AE18
Uniprot	P0AE18
GeneID	947882;
Calculated MW	30.4kDa
Target Sequence	AISIKTPEDIEKMRVAGRLAAEVLEMIEPYVKPGVSTGELDRICNDYIVNEQHAVSACLGYHGYPKSVCISINEV CHGIPDDAKLLKDGDIVNIDVTVIKDGFGDTSKMFIVGKPTIMGERLCRITQESLYLALRMVKPGINLREIGAAI QKFVEAEGFSVVREYCGHGIGRGFHEEPQVLHYDSRETNVVLKPGMTFTIEPMVNAGKKEIRTMKDGWTVKT KDRSLSAQYEHTIVWTDNGCEILTRKDDTIPAIISHDELEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 50% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only