

# Recombinant E. coli G,U Mismatch-Specific DNA Glycosylase,Mug (C-6His)

Catalog No: #AP74887

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74887-1 10ug #AP74887-2 50ug #AP74887-3 500ug #AP74887-4 1mg

## Description

Product Name	Recombinant E. coli G,U Mismatch-Specific DNA Glycosylase,Mug (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	E.coli
Accession No.	P0A9H1
Uniprot	P0A9H1
GeneID	947560;
Calculated MW	19.7kDa
Target Sequence	MVEDILAPGLRVVFCGINPGLSSAGTGFPFAHPANRFWKVIYQAGFTDRQLKPQEAQHLLDYRCGVTKLVDRP TVQANEVSKQELHAGGRKLIKIEDYQPQALAILGKQAYEQGFSQRGAQWGKQTLTIGSTQIWVLPNPSGLSR VSLEKLVEAYRELDQALVVRGRLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 2.5mM β-ME, 1mM PMSF, 50% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only