

Recombinant Human Brain Natriuretic Peptide, BNP (His27-His134, N-6His)



Catalog No: #AP74889

Orders: order@signalwayantibody.com

Package Size: #AP74889-1 10ug #AP74889-2 50ug #AP74889-3 500ug #AP74889-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Brain Natriuretic Peptide, BNP (His27-His134, N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16860
Uniprot	P16860
GeneID	4879;
Calculated MW	14.2kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMHLGSPGSASDLETSGLQEQRNHLQGKLSLQVEQTSLEPLQESPRPTGV WKSREVATEGIRGHRKMVLYTLRAPRSPKMQGSGCFGRKMDRISSSSGLGCKVLRH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 1mM DTT, 20% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only