

Recombinant Human Homeobox Protein OTX2 (C-6His)

Catalog No: #AP74890

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74890-1 10ug #AP74890-2 50ug #AP74890-3 500ug #AP74890-4 1mg

Description

Product Name	Recombinant Human Homeobox Protein OTX2 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P32243
Uniprot	P32243
GeneID	5015;
Calculated MW	33.47kDa
Target Sequence	MMSYLKQPPYAVNGLSLTTSGMDLLHPSVGYPGPWASCPAATPRKQRRERTTFTRAQLDVLEALFAKTRYPD IFMREEVALKINLPESRVQVWFKNRRRAKCRQQQQQQNGGQNKVVPKAKKTSPAREVSSESSTSGQFTPPS STSVPTIASSSAPVSIWSPASISPLSDPLSTSSSCMQRSPMTYTYQASGYSQGYAGSTSYFGGMDCGSYLTP MHHQLPGPGATLSPMGNTAVTSHLNQSPASLSTQGYGASSLGFNSTTDCLDYKDQTASWKLNFNADCLDYK DQTSSWKFQVLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only