

Recombinant Human Protein Phosphatase 1A, PPM1A (C-6His)

Catalog No: #AP74894

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74894-1 10ug #AP74894-2 50ug #AP74894-3 500ug #AP74894-4 1mg

Description

Product Name	Recombinant Human Protein Phosphatase 1A, PPM1A (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P35813
Uniprot	P35813
GeneID	5494;
Calculated MW	43.5kDa
Target Sequence	<p>MGAFLDKPKMEKHNAQQGGNGRLRYGLSSMQGWRVEMEDAHTAVIGLPSGLESWSFFAVYDGHAGSQVAKY</p> <p>CCEHLLDHITNNQDFKGSAGAPSVENVKNGIRTGFLEIDEHMRVMSEKKHGADRSGSTAVGVLISPQHTYFIN</p> <p>CGDSRGLLCRNKRVHFFTQDHKPSNPLEKERIQNAGGSVMIQRVNGSLAVSRALGDFDYKCVHGKGPTEQL</p> <p>VSPEPEVHDIERSEEDDQFIIACDGIWDVGMNEELCDFVRSRLEVTDDELEKVCNEVVDTCLYKGSRDNMSVIL</p> <p>ICFPNAPKVSPEAVKKEAELDKYLECRVEEIIKKQGEGVPDLVHVMRTLASENIPSLPPGGELASKRNVIEAVYN</p> <p>RLNPYKNDDTDSTSTDDMWLEHHHHHH</p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only