

Recombinant Human Protein Phosphatase 1 Regulatory Subunit 14A,PPP1R14A (C-6His)

Catalog No: #AP74895

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74895-1 10ug #AP74895-2 50ug #AP74895-3 500ug #AP74895-4 1mg

Description

Product Name	Recombinant Human Protein Phosphatase 1 Regulatory Subunit 14A,PPP1R14A (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96A00
Uniprot	Q96A00
GeneID	94274;
Calculated MW	17.76kDa
Target Sequence	MAAQRLGKRVLSKLGSPSRARGPGGSPGGLQKRHARVTVKYDRRELQRRLDVEKWIDGRLEELYRGMEAD MPDEINIDELLELESEEEERSRKIQGLLKSCGKPVEDFIQELLAKLQGLHRQPGLRQPSPSHDGSLSPLQDRART AHPLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 0.2mM EDTA, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only