

Recombinant Human Protein Phosphatase 1 Regulatory Subunit 1A,PPP1R1A (C-6His)

Catalog No: #AP74896

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74896-1 10ug #AP74896-2 50ug #AP74896-3 500ug #AP74896-4 1mg

Description

Product Name	Recombinant Human Protein Phosphatase 1 Regulatory Subunit 1A,PPP1R1A (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13522
Uniprot	Q13522
GeneID	5502;
Calculated MW	19.98kDa
Target Sequence	MEQDNSPQKIQFTVPLLEPHLDPEAAEQIRRRRPTPATLVLTSDQSSPEIDEDRIPNPHLKSTLAMSPRQRKKM TRITPTMKELQMMVEHHLGQQQQGEEPEGAAESTGTQESRPPGIPDTEVESRLGTSGTAKKTAECIPKThER GSKEPSTKEPSTHIPPLDSKGANSVLEHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 50% Glycerol, pH 8.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only