

## Recombinant Human Major Prion Protein, PRP, CD230

Catalog No: #AP74897



Package Size: #AP74897-1 10ug #AP74897-2 50ug #AP74897-3 500ug #AP74897-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Major Prion Protein, PRP, CD230
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04156
Uniprot	P04156
GeneID	5621;
Calculated MW	16.25kDa
Target Sequence	MGQGGGTHSQWNKPSKPKTNMKHMAGAAAAGAVVGGGLGGYVLGSAMSRPIIHFGSDYEDRYRENMHRY PNQVYYRPMDEYSNQNNFVHDCVNITIKQHTVTTTTKGENFTETDVKMMERVVEQMCITQYERESQAYYQRG SS
Formulation	Lyophilized from a 0.2 µm filtered solution of 5mM PB, 200mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only