

Recombinant Human Protein Kinase C Type ϵ , PKC ϵ , PKCE (C-6His)

Catalog No: #AP74898

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74898-1 10ug #AP74898-2 50ug #AP74898-3 500ug #AP74898-4 1mg

Description

Product Name	Recombinant Human Protein Kinase C Type ϵ , PKC ϵ , PKCE (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q02156
Uniprot	Q02156
GeneID	5581;
Calculated MW	19.64kDa
Target Sequence	MQELEYGPSVDWWALGVLMYEMMAGQPPFEADNEDDLFESILHDDVLYPVWLSKEAVSILKAFMTKNPHKRL GCVASQNGEDAIAKQHPFFKEIDWVLEQKKIKPPFKPRIKTKRDVNNFDQDFTRREEPVLTLDVAIVKQINQEEF KGFSYFGEDLMPLEHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only