

Recombinant Human Gankyrin,PSDM10

Catalog No: #AP74899



Package Size: #AP74899-1 10ug #AP74899-2 50ug #AP74899-3 500ug #AP74899-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Gankyrin,PSDM10
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75832
Uniprot	O75832
GeneID	5716;
Calculated MW	24.43kDa
Target Sequence	MEGCVSNLMVCNLAYSGLKLEELKESILADKSLATRTDQDSRTALHWACSAGHTEIVEFLLQLGVPVNDKDDAG WSPLHIAASAGRDEIVKALLGKGAQVNAVNQNGCTPLHYAASKNRHEIAVMLLEGGANPDAKDHYEATAMHR AAAKGNLKMIIHILLYKASTNIQDTEGNTPLHLACDEERVEEAKLLVSQGASIIENKKEETPLQVAKGGLGLILK RMVEG
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only