

# Recombinant Human Phosphoserine Phosphatase,PSP (C-6His)



Catalog No: #AP74900

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74900-1 10ug #AP74900-2 50ug #AP74900-3 500ug #AP74900-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Phosphoserine Phosphatase,PSP (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P78330
Uniprot	P78330
GeneID	5723;
Calculated MW	26.07kDa
Target Sequence	MVSHSELRKLFFYSADAVCFDVDSTVIREEGIDELAKICGVEDAVSEMTRRAMGGAVPFKAALTERLALIQPSRE QVQRLIAEQPPHLTPGIRELVSRQLQERNVQVFLISGGFRSIVEHVASKLNIPATNVFANRLKFYFNGEYAGFDET QPTAESGGKGVKIKLLKEKHFHFKIIMIGDGATDMEACPPADAFIFGGNVIRQQVKDNAKWYITDFVELLGELE ELEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 50mM Tris, 250mM NaCl, 50mM Imidazole, pH 8.5 .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only