

Recombinant Human 40S Ribosomal Protein S19,RPS19



Catalog No: #AP74903

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74903-1 10ug #AP74903-2 50ug #AP74903-3 500ug #AP74903-4 1mg

Description

Product Name	Recombinant Human 40S Ribosomal Protein S19,RPS19
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P39019
Uniprot	P39019
GeneID	6223;
Calculated MW	16.1kDa
Target Sequence	PGVTVKDVNQEFVRALAAFLKKSGLKVPEWVDTVKLAKHKELAPYDENWFYTRAASTARHLYLRGGAGV GSMTKIYGGRQRNGVMPSHFSRGSKSVARRVLQALEGLKMVEKDQDGGKRLTPQGQRDLDRAGQVAAAN KKH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only