

# Recombinant Human Sentrin-Specific Protease 7, SENP7 (N-6His)



Catalog No: #AP74904

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74904-1 10ug #AP74904-2 50ug #AP74904-3 500ug #AP74904-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Human Sentrin-Specific Protease 7, SENP7 (N-6His)  |
| Host Species    | E.coli   |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Human  |
| Accession No.   | Q9BQF6   |
| Uniprot         | Q9BQF6   |
| GeneID          | 57337;   |
| Calculated MW   | 27.3kDa  |
| Target Sequence | MKLKSVSQPSNTDAAKPTYTFLQKQSSGCYSLSITSNPDEEWREVRHTGLVQKLIVYPPPTKGGLGVTNEDL<br>ECLEEGEFLNDVIIDFYLYLILEKASDELVERSHIFSSFFYKCLTRKENNLTEDNPNL SMAQR RHKRVRTWTRH<br>INIFNKDYIFVPVNESSHWYLA |
| Formulation     | Supplied as a 0.2 µm filtered solution of 20mM HEPES, 5% Glycerol, pH 7.4.   |
| Storage         | Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.  |

Note: This product is for in vitro research use only