

Recombinant Human Sentrin-Specific Protease 8, SENP8 (N-6His)

Catalog No: #AP74905

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74905-1 10ug #AP74905-2 50ug #AP74905-3 500ug #AP74905-4 1mg

Description

Product Name	Recombinant Human Sentrin-Specific Protease 8, SENP8 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96LD8
Uniprot	Q96LD8
GeneID	123228;
Calculated MW	26.24kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMDPVVLSYMSLLRQSDVSLDPPSWLNDHIIGFAFEYFANSQFHDCSDHVS FISPEVTQFIKCTSNPAEIAMFLEPLDLPNKRVVFLAINDNSNQAAGGTHWLLVYLQDKNSFFHYDHSRSNS VHAKQVAEKLEAFLGRKGDKLAFVEEKAPAQQNSYDCGMYVICNTEALCQNFFRQQTESLLQLLTPAYITKKR GEWKDIATLAKK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only