

## Recombinant Human Stratifin,SFN

Catalog No: #AP74906



Package Size: #AP74906-1 10ug #AP74906-2 50ug #AP74906-3 500ug #AP74906-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Stratifin,SFN
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P31947
Uniprot	P31947
GeneID	2810;
Calculated MW	27.77kDa
Target Sequence	MERASLIQKAKLAEQAERYEDMAAFMKGAVEKGEELSCEERNLLSVAYKNVVGGQRAAWRVLSSIEQKSNEE GSEEEKGPEVREYREKVVETELQGVCDTVLGLLDSHLIKEAGDAESRVFYLMKMGDYRYRLAEVATGDDKKRIID SARSAYQEAMDISKKEMPPTNPIRLGLALNFSVFHYEIANSPPEAISLAKTTFDEAMADLHTLSEDSYKDSTLIM QLLRDNLTLWTADNAGEEGGEAPQEPQS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 250mM NaCl, 1mM EDTA, 1mM DTT, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only