

Recombinant Human Stathmin,STMN1 (C-6His)

Catalog No: #AP74907



Package Size: #AP74907-1 10ug #AP74907-2 50ug #AP74907-3 500ug #AP74907-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Stathmin,STMN1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16949
Uniprot	P16949
GeneID	3925;
Calculated MW	18.37kDa
Target Sequence	ASSDIQVKELEKRASGQAFELILSPRSKESVPEFPLSPPKKKDLSEEEIQKKLEAAEERRKSHEAEVLKQLAEKR EHEKEVLQKAIENNFSKMAEEKLTHKMEANKENREAQMAAKLERLREKDKHIEEVRKNKESKDPADETEA DLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only