

Recombinant Human Microtubule-Associated Protein Tau,MAPT-D,Tau-D (C-6His)



Catalog No: #AP74910

Orders: order@signalwayantibody.com

Package Size: #AP74910-1 10ug #AP74910-2 50ug #AP74910-3 500ug #AP74910-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Microtubule-Associated Protein Tau,MAPT-D,Tau-D (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P10636
Uniprot	P10636
GeneID	4137;
Calculated MW	15.37kDa
Target Sequence	MQIVYKPVDSLKVTSKCGSLGNIHHKPGGGQVEVKSEKLDKDRVQSKIGSLDNITHVPGGGNKKIETHKLTFR ENAKAKTDHGAEIVYKSPVSGDTSRHLNSVSTGSIDMVDSPQLATLADEVASLAKQLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only