

Recombinant Human Platelet-Derived Growth Factor BB,PDGF-BB



Catalog No: #AP74931

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74931-1 10ug #AP74931-2 50ug #AP74931-3 500ug #AP74931-4 1mg

Description

Product Name	Recombinant Human Platelet-Derived Growth Factor BB,PDGF-BB
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01127
Uniprot	P01127
GeneID	5155;
Calculated MW	12.42kDa
Target Sequence	MSLGSLTIAEPAMIAECKTRTEVFEISRRLIDRTNANFLVWPPCVEVQRCSGCCNNRNVQCRPTQVQLRPVQV RKIEIVRKKPIFKKATVTLEDHLACKCETVAAARPVT
Formulation	Lyophilized from a 0.2 µm filtered solution of 4mM HCl.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only