

Recombinant Human EIF4E-Binding Protein 1, EIF4EBP1 (N-6His)

Catalog No: #AP74932

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74932-1 10ug #AP74932-2 50ug #AP74932-3 500ug #AP74932-4 1mg

Description

Product Name	Recombinant Human EIF4E-Binding Protein 1, EIF4EBP1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13541
Uniprot	Q13541
GeneID	1978;
Calculated MW	14.74kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSGSSCSQTPSRAIPATRRVVLGDGVQLPPGDYSTTPGGTLFSTTPGGTR IYDRKFLMECRNSPVTKTPPRDLPTIPGVTSPSSDEPPMEASQSHLRNSPEDKRAGGEESQFEMDI
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only