

## Recombinant Human LMW-PTP,ACP1 (C-6His)

Catalog No: #AP74933



Package Size: #AP74933-1 10ug #AP74933-2 50ug #AP74933-3 500ug #AP74933-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human LMW-PTP,ACP1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P24666
Uniprot	P24666
GeneID	52;
Calculated MW	19.04kDa
Target Sequence	AEQATKSVL FVCLGNICRSPIAEAVFRKLVTDQNI SENWVIDSGAVSDWNVGRSPDPRAVSCLRNHGIHTAHK ARQITKEDFATFDYILCMDES NLRDLNRKSNQVK TCKAKIELLGSYDPQKQLIIEDPYYGNDSDFETVYQQCVR CCRAFLEKAHLEHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only