

Recombinant Human Annexin A13,ANXA13

Catalog No: #AP74936



Package Size: #AP74936-1 10ug #AP74936-2 50ug #AP74936-3 500ug #AP74936-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Annexin A13,ANXA13
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P27216
Uniprot	P27216
GeneID	312;
Calculated MW	35.4kDa
Target Sequence	GNRHAKASSPQGFVDVDRDAKLNKACKGMGTNEAAIIIEILSGRTSDERQQIKQKYKATYGKELEEVLKSEL SGNFEKTALALLDHPSEYAARQLQKAMKGLGTDESVLIEVLCTR TNKEIIAIKEAYQRLFDRSLES SDVKGDTSGNLK KILVSLLQANRNEGDDVDKDLAQDAKDLYDAGEGRWGTD ELAFNEVLAKRSYKQLRATFQAYQILIGKDIEE AIEEETSGDLQKAYLTLV RCAQDCEDYFAERLYKSMKGAGTDEETLIRIIVTRAEVDLQGIGAKAFQEKYQKSLSD MVRSDTSGDFRKL L VALLH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only