

## hnRNP K (Phospho-Ser284) Polyclonal Antibody

Catalog No: #12243

Package Size: #12243-1 50ul #12243-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	hnRNP K (Phospho-Ser284) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	Phospho-hnRNP K (S284) Polyclonal Antibody detects endogenous levels of hnRNP K protein only when phosphorylated at S284.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human hnRNP K around the phosphorylation site of S284.
Target Name	hnRNP K
Modification	Phospho
Other Names	HNRNPK; HNRPK; Heterogeneous nuclear ribonucleoprotein K; hnRNP K; Transformation up-regulated nuclear protein; TUNP
Accession No.	Swiss-Prot: P61978NCBI Gene ID: 3190
Uniprot	P61978
GeneID	3190;
Target Species	human
SDS-PAGE MW	51kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/20000

Not yet tested in other applications

## Images

Western Blot analysis of 3T3 cells using Phospho-hnRNP K (S284) Polyclonal Antibody



Note: This product is for in vitro research use only