

Recombinant Human Annexin A5, ANXA5

Catalog No: #AP74938



Package Size: #AP74938-1 10ug #AP74938-2 50ug #AP74938-3 500ug #AP74938-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Annexin A5, ANXA5
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08758
Uniprot	P08758
GeneID	308;
Calculated MW	35.9kDa
Target Sequence	MAQVLRGTVTDFPGFDERADAETLRKAMKGLGTDEESILTLTSSRSNAQRQEISAAFKTLFGRDLLDDLKSELT GKFEKLIVALMKPSRLYDAYELKHALKGAGTNEKVLTEIISRTPEELRAIKQVYEEYGSSEDDVVGDTSGYY QRMLVLLQANRDPDAGIDEAQVEQDAQALFQAGELKWTGDEEKFITIFGTRSVSHLRKVFVKYMTISGFQIEE TIDRETSNLEQLLAVVKSIRSIPAYLAETLYYAMKGAGTDDHTLIRVMVSRSETDLFNIRKEFRKNFATSLYS MIKGDTSGDYKALLLLCGEDD
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only