

Recombinant Human Annexin A6, ANXA6 (C-6His)

Catalog No: #AP74939



Package Size: #AP74939-1 10ug #AP74939-2 50ug #AP74939-3 500ug #AP74939-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Annexin A6, ANXA6 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08133
Uniprot	P08133
GeneID	309;
Calculated MW	76.9kDa
Target Sequence	AKPAQGAKYRGSIHDFPGFDPNQDAEALYTAMKGFSGDKEAILDIITSRSNRQRQEVCCQSYKSPLYGKDLIADLK YELTGKFERLIVGLMRPPAYCDAKEIKDAISGIGTDEKCLIEILASRTNEQMHQLVAAYKDAYERDLEADIIGDTS GHFQKMLVLLQGTREEDDVVSEDLVQQDVQDLYEAGELKWTDEAQFIYILGNRSKQHLRLVFDEYKTTG KPIEASIRGELSGDFEKLMLAVVKCIRSTPEYFAERLFKAMKGLGTRDNTLIRIMVSRSELDMLDIREIFRTKYEK SLYSMIKNDTSGEYKTKLLKLSGGDDDAAGQFFPEAAQVAYQMWELSAVARVELKGTVRPANDFNPDADAKA LRKAMKGLGTDEDTIIDITHRSNVQRQQIRQTFKSHFGRDLMTDLKSEISGDLARLILGLMMPPAHYDAKQLKK AMEGAGTDEKALIEILATRTNAEIRAINAEYKEDYHKSLEDALSSDTSGHFRRILISLATGHREEGGENLDQARE DAQVAEILEIADTPSGDKTSLETRFMILCTRSYPHLRRVFQEFIKMTNYDVEHTIKKEMSGDVRDAFVAIVQS VKNKPLFFADKLYKSMKGAGTDEKTLTRIMVSRSEIDLLNIRREFIEKYDKSLHQAIEGDTSGDFLKALLALCGG EDLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only