

Recombinant Human β -Ureidopropionase, UPB1 (C-6His)

Catalog No: #AP74942

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74942-1 10ug #AP74942-2 50ug #AP74942-3 500ug #AP74942-4 1mg

Description

Product Name	Recombinant Human β -Ureidopropionase, UPB1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UBR1
Uniprot	Q9UBR1
GeneID	51733;
Calculated MW	44.22kDa
Target Sequence	MAGAEWKSLEECLEKHLPLPDLQEVRVLYGKELRKLDLPREAFEAAASREDFELQGYAFEAAEEQLRRPRIVH VGLVQNRIPLPANAPVAEQVSALHRRIKAIVEVAAMCGVNIICFQEAWTMPFAFCTREKLPWTEFAESAEDGPT TRFCQKLAKNHDMVVVSPILERDSEHGDVLWNTAVVISNSGAVLGKTRKNHIPRVGDFNESTYYMEGNLGH VFQTQFGRIAVNICYGRHHPLNWL MYSINGAEIIFNPSATIGALSESLWPIEARNAAIANHCFTCAINRVGTEHFP NEFTSGDGKKAHQDFGYFYGSSYVAAPDSSRTPGLSRSDGLLVAKLDLNLCCQQVNDVWNFKMTGRYEMY ARELAEAVKSNYSPTIVKEVEHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only