

Recombinant Human Cyclin-Dependent Kinase 2,CDK2 (N-6His)

Catalog No: #AP74943

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74943-1 10ug #AP74943-2 50ug #AP74943-3 500ug #AP74943-4 1mg

Description

Product Name	Recombinant Human Cyclin-Dependent Kinase 2,CDK2 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P24941
Uniprot	P24941
GeneID	1017;
Calculated MW	36.1kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMENFQKVEKIGEGTYGVVYKARNKLTGEVVALKKIRLDTETEGVPSTAIREIS LLKELNHPNIVKLLDVIHTENKLYLVFEFLHQDLKKFMDASALTGIPLPLIKSYLFQLLQGLAFCHSHRVLHRDLK PQNLLINTEGAIKLADFGLARAFGVPVRTYTHEVVTWYRAPEILLGCKYYSTAVDIWSLGCIFAEMVTRRALFP GDSEIDQLFRIFRTLGTPEVVWPGVTSMPDYKPSFPKWARQDFSKVVPPLDEDGRSLLSQMLHYDPNKRIS AKAALAHPPFQDVTKPVPHLRL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, 40% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only