

Recombinant Human C-Reactive Protein,CRP (C-6His)

Catalog No: #AP74945

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74945-1 10ug #AP74945-2 50ug #AP74945-3 500ug #AP74945-4 1mg

Description

Product Name	Recombinant Human C-Reactive Protein,CRP (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02741
Uniprot	P02741
GeneID	1401;
Calculated MW	23.04kDa
Target Sequence	MQTDMSRKAFVFPKESDTSYVSLKAPLTKPLKAFVCLHFDYELSSSTRGYSIFSYATKRQDNEILIFWSKDIGYS FTVGGSEILFEVPEVTVAPVHICTSWESASGIVEFWVDGKPRVRKSLKKGTVGAEASIILGQEQDSFGGNFEG SQSLVGDIGNVMWDFVLSPEINTIYLGPFSPNVLNWRALKYEVQGEVFTKPQLWPLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mMTris2MUreapH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only