

Recombinant Human γ-Glutamylcyclotransferase,GGCT (N-6His)

Catalog No: #AP74946

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74946-1 10ug #AP74946-2 50ug #AP74946-3 500ug #AP74946-4 1mg

Description

Product Name	Recombinant Human γ-Glutamylcyclotransferase,GGCT (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75223
Uniprot	O75223
GeneID	79017;
Calculated MW	23.17kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMANS GCKDVTGPDEESFLYFAYGSNLLTERIHLRNPSAAFFCVARLQDFKLD FGNSQGKTSQTVHGGIATIFQSPGDEVWGVVWKMNKS NLNSLDEQEGVKSGMYV VIEVKVATQEGKEITCR SYLMTNYESAPPSPQYKKIICMGAKENGLPLEYQEKLKAIEPNDYTGKVSEEIEDI IKKGETQTL
Formulation	Supplied as a 0.2 μm filtered solution of 20mM TrisHCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only