

Recombinant Human Cyclophilin C, PPlase C, PPIC (N-Trx, 6His)



Catalog No: #AP74947

Orders: order@signalwayantibody.com

Package Size: #AP74947-1 10ug #AP74947-2 50ug #AP74947-3 500ug #AP74947-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cyclophilin C, PPlase C, PPIC (N-Trx, 6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P45877
Uniprot	P45877
GeneID	5480;
Calculated MW	33.78kDa
Target Sequence	MSDKIIHLTDDSFDTDLKADGAILVDFWAEWCGPCKMIAPILDEIADEYQGKLTVAKLNIDQNPGTAPKYGIRGI PTLLLFKNGEVAATKVGALSKGQLKEFLDANLAGSGSGMHMHHHHSSGLVPRGSGMKETA AAKFERQHMD SPDLGTDDDDKAMAKRGPSVTAKVFFDVRIGDKDVGRIVIGLFGKVVPKTVENFVALATGEKGYGYKGSKFHR VIKDFMIQGGDITGDGTGGVSIYGETFPDENFKLKHYGIGWVSMANAGPDTNGSQFFITLTKPTWLDGKHVV FGKVIDGMTVVHSIELQATD
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only