

# Recombinant Human Cyclophilin E, PPIase E, PPIE (N-6His)

Catalog No: #AP74948

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74948-1 10ug #AP74948-2 50ug #AP74948-3 500ug #AP74948-4 1mg

## Description

Product Name	Recombinant Human Cyclophilin E, PPIase E, PPIE (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UNP9
Uniprot	Q9UNP9
GeneID	10450;
Calculated MW	35.6kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMATTKRVLVVGGLAEEVDDKVLHAAFIPFGDITDIQIPLDYETEKHRGF AFVE FELAEDAAAAIDNMNESELFGRITIRVNLAKPMRIKEGSSRPVWSDDDWLKKFSGKTLEENKEEEGSEPPKAET QEGEPIAKKARSNPQVYMDIKIGNKPAGRIQMLLRSDVVPMTAENFRCLCTHEKGF GFKGSSFHRIIPQFMCQ GGDFTNHNGTGGKSIYGKKFDDENFILKHTGPGLLSMANSGPNTNGSQFFLTCDKTDWLDGKHVVFGEVTEG LDVLRQIEAQGSKDGKPKQKVIADCGEYV
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only