

## Recombinant Human Cytochrome b5,CYB5A (N-6His)

Catalog No: #AP74949



Package Size: #AP74949-1 10ug #AP74949-2 50ug #AP74949-3 500ug #AP74949-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Cytochrome b5,CYB5A (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00167
Uniprot	P00167
GeneID	1528;
Calculated MW	17.5kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAEQSDEAVKYYTLEEIQKHNSKSTWLILHHKVYDLTKFLEEHPGGEEVLR EQAGGDATENFEDVGHSTDAREMSKTFIIGELHPDDRPKLNKPPETLITTIDSSSSWWTNWWIPAISAVAVALM YRLYMAED
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 0.1mM EDTA, pH 7.25.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only