

Recombinant Human Cytoplasmic Dynein Light Chain 1,DYNLL1 (N-6His)

Catalog No: #AP74950

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74950-1 10ug #AP74950-2 50ug #AP74950-3 500ug #AP74950-4 1mg

Description

Product Name	Recombinant Human Cytoplasmic Dynein Light Chain 1,DYNLL1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P63167
Uniprot	P63167
GeneID	8655;
Calculated MW	12.5kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMCDRKAIVIKNADMSEEMQQDSVECATQALEKYNIKDIAAHIKKEFDKKYNP TWHCIVGRNFGSYVTHEKHFYFYLGQVAILLFKSG
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only