

Recombinant Human EIF1A, X-Chromosomal, EIF1AX (N-6His)



Catalog No: #AP74951

Orders: order@signalwayantibody.com

Package Size: #AP74951-1 10ug #AP74951-2 50ug #AP74951-3 500ug #AP74951-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EIF1A, X-Chromosomal, EIF1AX (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P47813
Uniprot	P47813
GeneID	1964;
Calculated MW	18.6kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMPKNGKGGKNNRRRGKNENESEKRELVFKEDGQEYAQVIKMLGNRLEAM CFDGVKRLCHIRGKLRKKVWINTSDIILVGLRDYQDNKADVILKYNADERSLKAYGELPEHAKINETDTFGPGD DDEIQFDDIGDDEDEDIDI
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only