

Recombinant Human Fatty Acid-Binding Protein 6, FABP6, I-BABP (N-6His)

Catalog No: #AP74953

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74953-1 10ug #AP74953-2 50ug #AP74953-3 500ug #AP74953-4 1mg

Description

Product Name	Recombinant Human Fatty Acid-Binding Protein 6, FABP6, I-BABP (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51161
Uniprot	P51161
GeneID	2172;
Calculated MW	16.6kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAFTGKFEMESEKNYDEFMKLLGISSDVIEKAHNFKIVTEVQQDGDFTWS QHYYGGHTMTNKFTVGKESNIQTMGGKTFKATVQMEGGKLVNFPNYHQTSEIVGDKLVEVSTIGGVTYERV SKRLA
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 0.5mM DTT, 50% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only