

# Recombinant Human Fibroblast Growth Factor 12,FGF-12



Catalog No: #AP74954

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74954-1 10ug #AP74954-2 50ug #AP74954-3 500ug #AP74954-4 1mg

## Description

Product Name	Recombinant Human Fibroblast Growth Factor 12,FGF-12
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P61328
Uniprot	P61328
GeneID	2257;
Calculated MW	20.45kDa
Target Sequence	MESKEPQLKGIVTRLFSQQGYFLQMHPDGTIDGTDKDNSDYTLFNLIPVGLRVVAIQGVKASLYVAMNGEGYL YSSDVFTPECKFKESVFENYYVIYSSTLYRQQESGRAWFLGLNKEGQIMKGNRVKKTTPSSHFPKPIEVCMY REQSLHEIGEKQGRSRKSSGTPTMNGGKVVNQDST
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5mM EDTA, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only