

# Recombinant Human Fibroblast Growth Factor 21,FGF-21 (N-6His)

Catalog No: #AP74955

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74955-1 10ug #AP74955-2 50ug #AP74955-3 500ug #AP74955-4 1mg

## Description

Product Name	Recombinant Human Fibroblast Growth Factor 21,FGF-21 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NSA1
Uniprot	Q9NSA1
GeneID	26291;
Calculated MW	21.69kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMHPIPDSSPLLQFGGQVRQRYLYTDDAQQTEAHLEIREDGTVGGAADQSPE SLLQLKALKPGVIQILGVKTSRFLCQRPDGALYGSLHFDPEACSFRELLLEDGYNVYQSEAHGLPLHLPGNKSP HRDPAPRGPAPRFLPLPGLPPAPPEPPGILAPQPPDVGSSDPLSMVGPSQGRSPSYAS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, 2mM EDTA, pH9.0 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only