

Recombinant Human Fumarase, FH

Catalog No: #AP74956



Package Size: #AP74956-1 10ug #AP74956-2 50ug #AP74956-3 500ug #AP74956-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fumarase, FH
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07954
Uniprot	P07954
GeneID	2271;
Calculated MW	50.2kDa
Target Sequence	MASQNSFRIEYDTFGELKVPNDKYYGAQTVRSTMFKIGGVTERMPTPVKAFGILKRAAAEVNQDYGLDPKIA NAIMKAADEVAEGKLNDFPLVWVWQTGSGTQTNMNVNEVISNRAIEMLGELGSKIPVHPNDHVNKSQSSND TFPTAMHIAAAIEVHEVLLPGLQKLHDALDAKSKEFAQIIKIGRTHQDAVPLTLGQEFSGYVQQVKYAMTRIKA AMPRIYELAAGGTAVGTGLNTRIGFAEKVAAKVAALTGLPFVTAPNKFEALAAHDALVELSGAMNTTACSLMKI ANDIRFLGSGPRSGLGELILPENEPGSSIMPGKVNPTQCEAMTMVAAQVMGNHVAVTVGGSNHGFELNVFKP MMIKNVLHSARLLGDASVSFTENCVVGIQANTERINKLMNESLMLVTALNPHIGYDKAAKIAKTAHKNGSTLKET AIELGYLTAEQFDEWVKPKDMLGPK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only